

MARTCARE® FACESHIELD SPARE VISOR

PRODUCT CODE: ANM060-730-000

REPLACEMENT POLYCARBONATE VISOR FOR THE MARTCARE® FACESHEILD AND THE PREMIUM MARTCARE® FACESHIELD

TECHNICAL
SPECIFICATION

MARTCARE® RANGE



INSTRUCTIONS FOR USE

- Before use, the browguard should be fitted to the visor
- In use, ensure that the faceshield fits comfortably and securely on the head
- Ensure that the visor is in place before entering the hazard area

MATERIALS

Visor: Polycarbonate

DIMENSIONS OF OUTER PACKAGING

- Quantity: 100
- Height: 140mm
- Width: 235mm
- Length: 420mm
- Weight: 9.48Kg (Approx.)

FACESHIELD FEATURES & PERFORMANCE

- Offers protection against medium energy high speed particles at 120m/s (270mph), liquid droplets and molten metal
- The optical qualities of the visor comply with the requirements of optical class 1 (highest)

CLEANING & MAINTENANCE

- The browguard may be cleaned with soap and warm water and dried with a soft cloth
- Do not clean with chemical or abrasive cleaners
- The eye shield should be stored out of direct sunlight, away from chemicals and not exposed to extremes of temperature

CONFORMITY

- Meets the requirements of EN 166.1.B.3.9 and EN 169
- EC Certificate Number 2406, issued by BSi Product Certification

DIMENSIONS OF VISOR

- Height: 210mm
- Width: 384mm
- Length: 1mm
- Weight: 82g (Approx.)

LIMITATIONS OF USE

- Protection will only be offered if the faceshield is fitted correctly
- If the faceshield becomes scratched or damaged they should be replaced
- Do not apply paint, solvents, adhesives or self-adhesive labels, except in accordance with instructions from the manufacturer
- The faceshield, when in contact with the skin may cause allergic reactions to susceptible individuals. If this occurs leave the hazard area, remove the faceshield and seek medical advice
- Polycarbonate visors are not suitable for use with all chemicals - if you are unsure contact your supplier or JSP for help and advice

